

DIASYS 3 Plus



THE MOST ADVANCED HYBRID ABPM DELIVERING SUPERIOR ACCURACY

TRIPLE-SENSOR TECHNOLOGY (3 sensors: Microphone, Accelerometer and Barometer) delivering:

- Hybrid Auscultatory & Oscillometric Measurement
- Activity and Position Assessment
- Orthostatic Hypotension Assessment
- HR irregularity detection



Dual Sensor Wireless Patch adding these

UNIQUE FEATURES

- 24hr Arterial Compliance / Qkd
- Central Blood Pressure
- Enhanced Position Information
- Atrial Fibrillation Detection

Fully compatible with
our advanced software

HOLTERSOFTE

ULTIMA

Diasys 3 Plus

SPECIFICATIONS	Size	85x85x25mm
	Weight	161g
	Weight (inc. batteries)	195g
	IP Rating	IP22
	Mode	Hybrid Auscultatory + Oscillometric
	Position Sensor	2 sensors (an accelerometer and barometer)
	Cuff Connection	Instant connection and disconnection
	PC Connection	USB Type-C
	AF Feature	Screening

Optional accessories

Pouch and belt or belt dip, Range of Cuff Sizes, Standard or rechargeable batteries
Standard USB interface cable (USB Type-C)

Dual Sensor Wireless Patch

SPECIFICATIONS	Size	55x40x17mm
	Weight	23g
	Weight (inc. batteries)	52g
	IP Rating	IP68
	Sensors	ECG & Accelerometer
	AF Feature	ECG Based Detection

Önskar ni att vi kommer ut och presenterar produkten eller skickar en offert, kontakta vänligen:

NorDiaTech er eksklusiv forhandler i
Danmark og Sverige for:



4, passage Saint Antoine - 92500 Rueil Malmaison - France
Tél. +33 (0)1 41 39 01 61 - Fax +33 (0)1 41 39 01 99

NOVACOR is a registered ISO 9001 and ISO 13485 company.

Ulf Ilar

Produktspecialist / Product specialist

+46 70 366 88 13

ulf.ilar@nordiatech.dk

NorDiaTech

