

DIASYS 3



**THE ULTIMATE OSCILLOMETRIC ABPM
WITH VALIDATED UNPARALLELED ACCURACY**

A **Revolution** in Precision and Comfort

Oscillometric Measurement

+

Accelerometer for
Activity Assessment

□

Orthostatic Hypotension
Assessment

□

HR Irregularity Detection

Fully compatible with
our advanced software

HOLTERSOFTE

ULTIMA

Diasys 3

SPECIFICATIONS	Size	85x85x25mm
	Weight	158g
	Weight (inc. batteries)	192g
	IP Rating	IP22
	Mode	Oscillometric
	Position Sensor	Embedded multi axis accelerometer
	Cuff Connection	Instant connection and disconnection
	PC Connection	USB Type-C
	Arrhythmia Detection	Heart rate irregularity

Accessories include

Pouch and belt or belt clip, Range of Cuff Sizes, Standard or rechargeable batteries
Standard USB interface cable (USB Type-C).

Önskar ni att vi kommer ut och presenterar produkten eller skickar en offert, kontakta vänligen:

NorDiaTech er eksklusiv forhandler i
Danmark og Sverige for:


NOVACOR

4, passage Saint Antoine - 92500 Rueil Malmaison - France
Tél. +33 (0)1 41 39 01 61 - Fax +33 (0)1 41 39 01 99

NOVACOR is a registered ISO 9001 and ISO 13485 company.

Ulf Ilar

Produktspecialist / Product specialist

+46 70 366 88 13

ulf.ilar@nordiatech.dk

NorDiaTech

