

Daavlin

Let's Be Clear

Daavlin Phototherapy Solutions

DermaPal



The Daavlin DermaPal

The Daavlin DermaPal is ideal for scalp and spot treatments.
DermaPal is safe, affordable and powerful.

The Daavlin DermaPal IS SAFE, AFFORDABLE AND POWERFUL

The Daavlin DermaPal delivers a clear solution for psoriasis, eczema and vitiligo.

Effective

- Treatments last for a few seconds to a few minutes, about three times a week.
- Psoriasis patients may see a reduction in symptoms such as redness and flaking in as little as 6-8 treatments.
- Vitiligo patients may see repigmentation begin within 8 weeks.

Affordable

- Daavlin products are reasonably priced.
- In-office phototherapy treatments are reimbursable by Medicare and most insurance.

Safe and easy to use

- Safe for most patients, including pregnant women and children, as well as elderly and immuno-compromised patients.
- Plugs into a standard electrical outlet.
- DermaPal has a self-contained digital timer and doesn't require an inconvenient separate timer.



Patient Friendly

The astonishingly light-weight treatment wand provides treatment without tiring the user. The durable slide-on comb attachment also doubles as a distance guide. UV protective eyewear and a convenient carrying case are included.



Daavlin DermaPal Data

Treatment Area	1" x 4.5"	(2.54 cm x 11.4 cm)
Wand Size	8.25" x 1.75"	(21 cm x 5 cm)
Wand Weight	4 oz.	(0.1 kg)
Fixture Size	9" x 5"	(23 cm x 13 cm)
Total Height	2.75"	(7 cm)
Total Weight:	1 lb 12 oz.	(0.8 kg)
Available in 110-120 V, 60 Hz and 220-240 V, 50 or 60 Hz.		

For further information please contact Daavlin:

Phone: 1.800.322.8546 Fax: 419.636.1739

International: 001.419.636.6304 Website: www.daavlin.com

All Daavlin phototherapy units are 510k cleared by the FDA.