

# Daavlin

## Let's Be Clear

Daavlin *Phototherapy Solutions*

### UVseries



### The Daavlin UV Series

The Daavlin UV Series is powerful and versatile, specifically designed for offices with limited space.

**The Daavlin UV Series IS ACCURATE , AFFORDABLE, VERSATILE AND COMPACT**

The Daavlin UV Series combines space-saving design with maximum efficiency. Twenty-four high output lamps in a compact space means short treatment times.

**Accurate**

Daavlin's integrating dosimetry system allows the operator to enter the dose in Joules for UVA and/or millijoules for UVB. The built-in metering system then automatically adjusts the treatment time to compensate for variations in output, delivering the exact dose every time. (Digital timer systems are available upon request).

Additional features include a digital lock out to prevent unauthorized use and the display of real-time power output readings on demand.

UV Series devices are manufactured with the control system located outside of the unit to prevent the patient from tampering with the control system while the treatment is in process.

**Affordable**

- The 120V (220V overseas) electrical requirements of the UV Series make it economical to operate, costing just pennies per treatment.
- Relamping the UV Series is less expensive as well, with only half the number of lamps of larger full-body cabinets.

**Versatile**

- The UV Series can be equipped with Narrow Band UVB, Broad Band UVB, or UVA lamps. Combination units are also available.
- A platform is available to elevate the patient to the height of the lamps for optimum treatment of the lower legs.
- Optional cooling fans directed downward from the top of the device make treatments more comfortable.

**Compact**

- The unique, elliptical design of our UV Series enhances the uniform irradiance of the device and also makes efficient use of valuable floor space in offices where space is a concern.
- With its doors folded inward, the UV Series can be easily rolled through a 27 inch (69cm) doorway.



**Daavlin UV Series Data**

Height	79 in	(195.6 cm)
Width	41 in	(104.1 cm)
Depth	30.5 in	(77.5 cm)
Depth (with doors open)	42 in	(106.7 cm)
Weight (24 lamps)	327 lbs	(148.3 kg)
Available in 110-120 V, 60 Hz and 220-240 V, 50 or 60 Hz.		

**For further information please contact Daavlin:**

Phone: 1.800.322.8546 Fax: 419.636.1739

International: 001.419.636.6304 Website: [www.daavlin.com](http://www.daavlin.com)

All Daavlin phototherapy units are 510k cleared by the FDA.