



EMODIAL



EMODIAL

Via Finati, 4/A - Ferrara - ITALY

Tel. +39 0532 978198 - Fax +39 0532 978414

www.emodial.com - info@emodial.com



Ref.GB003 - Rev. 07/2019

EMOSTAN[®]

Fistula Pad



EMOSTAN[®]

Fistula Pad

DESCRIPTION

EMOSTAN is a highly absorbent pad with compressive hemostatic action, specific for post-dialysis arteriovenous access care to be used after the removal of the fistula needles.

Also suitable for patients who require intravenous treatment and blood transfusion.

INTENDED USE

Post dialysis dressing to stop the bleeding after the removal of the fistula needles.

PURPOSE

Designed to stop bleeding and to reduce the haemorrhagic risk, thanks to the compression on the needle holes and to work as a barrier against the risk of infection.

The non-adherent, biocompatible layer on the pad doesn't stick to the skin and **guarantees a non-traumatic removal**.



DIRECTIONS FOR USE





REFERENCES

Data and technical specifications can be changed without notice.

CODE	PAD (cm)		PCS / PACK	PACKS / BOX
	Dimension	Thickness		
S0210	ø 2.6	0.8 ± 10%	2	500
S2010-1	ø 2.6	0.8 ± 10%	1	700
S0211	ø 3.5	0.8 ± 10%	2	500
S2011-1	ø 3.5	0.8 ± 10%	1	700
S0212	3.5 x 4.5	0.8 ± 10%	2	300
S2012-1	3.5 x 4.5	0.8 ± 10%	1	500

CLASSIFICATION • CLASS I MEDICAL DEVICE
BIOCOMPATIBLE • ACCORDING TO UNI EN ISO 10993 STANDARDS
STERILIZATION • GAMMA RAYS
EXPIRATION • 5 YEARS

  • SINGLE-USE